

## SO-LOW MODEL CH25-25

### CHEST STYLE FREEZER TO -25°C

<b>Capacity [ Cu. Ft. ]</b>	<b>25 Cubic Feet</b>
<b>Capacity [ Liters ]</b>	<b>707 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-0°C to -25°C</b>
<b>Temperature Range [ °F ]</b>	<b>+32°F to -13°F</b>

Dimensions (W x D x H)

Exterior [ Inches ] 83.25" x 32" x 33.75"

Exterior [ Centimeters ] 211 x 81 x 85

Interior [ Inches ] 76.75" x 23.5" x 24.88

Interior [ Centimeters ] 194 x 59 x 63

#### Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage 15 Amp Dedicated Line

2 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

#### Construction:

Temperature Control Digital Display  
High / Low Alarm  
Min / Max Temperature

Alarm Relay Dry Contacts

Defrost Manual

Opening Top Opening Lid

Lockable Lid/Door No

Interior / Exterior Steel

Insulation Polyurethane Foam

Compressors Hermetic Compressor

Primary Refrigerant(s) R134A

Air Circulation Gravity Flow


Casters 4, Lockable

Shelving None

Approximate BTU/Hour 600

Approximate Decibels 55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 1 Year Parts & 90 Days Labor