

## SO-LOW MODEL CH25-3

### MINI CHEST STYLE FREEZER TO -25°C

<b>Capacity [ Cu. Ft. ]</b>	<b>3.5 Cubic Feet</b>
<b>Capacity [ Liters ]</b>	<b>99 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-0°C to -25°C</b>
<b>Temperature Range [ °F ]</b>	<b>+32°F to -13°F</b>

Dimensions (W x D x H)

Exterior [ Inches ]	22" x 23.5" x 34"
Exterior [ Centimeters ]	56 x 59 x 86
Interior [ Inches ]	15.75" x 17" x 25.87"
Interior [ Centimeters ]	40 x 43 x 65

#### Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage	15 Amp Dedicated Line 2 Running Amps
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

#### Construction:

Temperature Control	Digital Display High / Low Alarm Min / Max Temperature
Alarm Relay	Dry Contacts
Defrost	Manual
Opening	Top Opening Lid
Lockable Lid/Door	No
Interior / Exterior	Steel
Insulation	Polyurethane Foam
Compressors	Hermetic Compressor
Primary Refrigerant(s)	R600
Air Circulation	Gravity Flow
Casters	4, Lockable
Shelving	None
Approximate BTU/Hour	600
Approximate Decibels	55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 1 Year Parts & 90 Days Labor