

## SO-LOW MODEL CH25-3

### MINI CHEST STYLE FREEZER TO -25°C

#### Specifications:

**Capacity** 3.5 Cu. Ft. | 99 Liters

**Temperature Range °C** -0°C to -25°C

**Temperature Range °F** +32°F to -13°F

Dimensions (W x D x H)

Exterior [ Inches ] 22" x 23.5" x 34"

Exterior [ Centimeters ] 56 x 59 x 86

Interior [ Inches ] 15.75" x 17" x 25.87"

Interior [ Centimeters ] 40 x 43 x 65

#### Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage 15 Amp Dedicated Line

2 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

#### Construction:

Temperature Control Digital Display  
High / Low Alarm  
Min / Max Temperature

Alarm Relay Dry Contacts

Defrost Manual

Opening Top Opening Lid

Lockable Lid/Door No

Interior / Exterior Steel

Insulation Polyurethane Foam

Compressors Hermetic Compressor

Primary Refrigerant(s) R600

Air Circulation Gravity Flow

Casters 4, Lockable

Shelving None

Approximate BTU/Hour 600

Approximate Decibels 55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



#### Warranty:

1 Year Parts / 90 Days Labor