

SO-LOW MODEL CH25-5

MINI CHEST STYLE FREEZER TO -25°C

Specifications:

Capacity 5 Cu. Ft. | 142 Liters

Temperature Range °C -0°C to -25°C

Temperature Range °F +32°F to -13°F

Dimensions (W x D x H)

Exterior [Inches] 29" x 24" x 33.63"

Exterior [Centimeters] 73 x 60 x 85

Interior [Inches] 24" x 15.25" x 26"

Interior [Centimeters] 60 x 38 x 66

Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage 15 Amp Dedicated Line

2 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

Construction:

Temperature Control Digital Display
High / Low Alarm
Min / Max Temperature

Alarm Relay Dry Contacts

Defrost Manual

Opening Top Opening Lid

Lockable Lid/Door No

Interior / Exterior Steel

Insulation Polyurethane Foam

Compressors Hermetic Compressor

Primary Refrigerant(s) R600

Air Circulation Gravity Flow

Casters 4, Lockable

Shelving None

Approximate BTU/Hour 600

Approximate Decibels 55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty:

1 Year Parts / 90 Days Labor