

SO-LOW MODEL CH25-9

CHEST STYLE FREEZER TO -25°C

Capacity [Cu. Ft.]	9 Cubic Feet
Capacity [Liters]	255 Liters
Temperature Range [°C]	-0°C to -25°C
Temperature Range [°F]	+32°F to -13°F

Dimensions (W x D x H)	
Exterior [Inches]	39" x 23.5" x 33.5"
Exterior [Centimeters]	99 x 59 x 85
Interior [Inches]	33.5" x 17.5" x 27.5"
Interior [Centimeters]	85 x 44 x 69

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage	15 Amp Dedicated Line 2 Running Amps
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

Construction:

Temperature Control	Digital Display High / Low Alarm Min / Max Temperature
Alarm Relay	Dry Contacts
Defrost	Manual
Opening	Top Opening Lid
Lockable Lid/Door	No
Interior / Exterior	Steel
Insulation	Polyurethane Foam
Compressors	Hermetic Compressor
Primary Refrigerant(s)	R600
Air Circulation	Gravity Flow
Casters	4, Lockable
Shelving	None
Approximate BTU/Hour	600
Approximate Decibels	55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & 90 Days Labor