

SO-LOW MODEL CH40-15

CHEST STYLE FREEZER TO -40°C

Capacity [Cu. Ft.]	15 Cubic Feet
Capacity [Liters]	424 Liters
Temperature Range [°C]	-0°C to -40°C
Temperature Range [°F]	+32°F to -40°F

Dimensions (W x D x H)	
Exterior [Inches]	55.5" x 29.75" x 34.38"
Exterior [Centimeters]	140 x 75 x 87
Interior [Inches]	49.25" x 21.63" x 25.75"
Interior [Centimeters]	125 x 59 x 65

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage	15 Amp Dedicated Line 3 Running Amps
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

Construction:

Temperature Control	Digital Display High / Low Alarm Min / Max Temperature
Alarm Relay	Dry Contacts
Defrost	Manual
Opening	Top Opening Lid
Lockable Lid/Door	No
Interior / Exterior	Steel
Insulation	Polyurethane Foam
Compressors	Hermetic Compressor
Primary Refrigerant(s)	R404A
Air Circulation	Gravity Flow
Casters	4, Lockable
Shelving	None
Approximate BTU/Hour	700
Approximate Decibels	55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & 90 Days Labor