

## SO-LOW MODEL CH40-25

### CHEST STYLE FREEZER TO -40°C

#### Specifications:

**Capacity** 25 Cu. Ft. | 707 Liters

**Temperature Range °C** -0°C to -40°C

**Temperature Range °F** +32°F to -40°F

Dimensions (W x D x H)

Exterior [ Inches ] 83.25" x 32" x 33.75"

Exterior [ Centimeters ] 211 x 81 x 85"

Interior [ Inches ] 76.75" x 23.5" x 24.88"

Interior [ Centimeters ] 194 x 59 x 63

#### Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage 15 Amp Dedicated Line

3 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

#### Construction:

Temperature Control Digital Display  
High / Low Alarm  
Min / Max Temperature

Alarm Relay Dry Contacts

Defrost Manual

Opening Top Opening Lid

Lockable Lid/Door No

Interior / Exterior Steel

Insulation Polyurethane Foam

Compressors Hermetic Compressor

Primary Refrigerant(s) R404A

Air Circulation Gravity Flow

Casters 4, Lockable

Shelving None

Approximate BTU/Hour 700

Approximate Decibels 55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



#### Warranty:

1 Year Parts / 90 Days Labor