

## *SO-LOW CHEST STYLE FREEZER TO -40°C*

### *MODEL CH40-25*

**25 Cubic Ft. (707 liters) Chest – Temperature Range: 0°C to -40°C**

#### **Top Opening**

#### ***Specifications:***

* Capacity	25 cubic ft.
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior w x d x h	83.25" x 32" x 33.75"
* Interior w x d x h	76.75" x 23.5" x 24.88"

#### **CONTROL**

Digital Temperature control with two Temperature displays. One display shows the set temperature, the other shows the actual temperature in the freezer.

#### **ALARM**

Battery operated alarm system provides An audible and visual signal in a high/low temperature alarm condition. The alarm will also signal in the event of a power failure. Also provided is an alarm silencing and battery test feature.

MANUAL DEFROST

CASTERS FOR EASY MOBILITY

HERMETIC MOTOR COMPRESSOR

115 VOLT, 60 HERTZ, 1 PHASE

WARRANTY One Year Parts  
90 Days Labor – Shipping Weight 450 lbs.

#### **TECHNICAL SUPPORT**

Supplied by Nation-Wide-Network of service companies.

