

SO-LOW MODEL CH40-3

MINI CHEST STYLE FREEZER TO -40°C

Specifications:

Capacity 3.5 Cu. Ft. | 99 Liters

Temperature Range °C -0°C to -40°C

Temperature Range °F +32°F to -40°F

Dimensions (W x D x H)

Exterior [Inches] 22" x 23.5" x 34"

Exterior [Centimeters] 55 x 59 x 86

Interior [Inches] 15.75" x 17" x 25.87"

Interior [Centimeters] 40 x 43 x 65

Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage 15 Amp Dedicated Line

2 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

Construction:

Temperature Control Digital Display
High / Low Alarm
Min / Max Temperature

Alarm Relay Dry Contacts

Defrost Manual

Opening Top Opening Lid

Lockable Lid/Door No

Interior / Exterior Steel

Insulation Polyurethane Foam

Compressors Hermetic Compressor

Primary Refrigerant(s) R290

Air Circulation Gravity Flow

Casters 4, Lockable

Shelving None

Approximate BTU/Hour 600

Approximate Decibels 55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty:

1 Year Parts / 90 Days Labor