

SO-LOW MODEL CH40-3

MINI CHEST STYLE FREEZER TO -40°C

Capacity [Cu. Ft.]	3.5 Cubic Feet
Capacity [Liters]	99 Liters
Temperature Range [°C]	-0°C to -40°C
Temperature Range [°F]	+32°F to -40°F

Dimensions (W x D x H)

Exterior [Inches]	22" x 23.5" x 34"
Exterior [Centimeters]	55 x 59 x 86
Interior [Inches]	15.75" x 17" x 25.87"
Interior [Centimeters]	40 x 43 x 65

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage	15 Amp Dedicated Line 2 Running Amps
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

Construction:

Temperature Control	Digital Display High / Low Alarm Min / Max Temperature
Alarm Relay	Dry Contacts
Defrost	Manual
Opening	Top Opening Lid
Lockable Lid/Door	No
Interior / Exterior	Steel
Insulation	Polyurethane Foam
Compressors	Hermetic Compressor
Primary Refrigerant(s)	R290
Air Circulation	Gravity Flow
Casters	4, Lockable
Shelving	None
Approximate BTU/Hour	600
Approximate Decibels	55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & 90 Days Labor