

## SO-LOW MODEL CH40-9

### CHEST STYLE FREEZER TO -40°C

<b>Capacity [ Cu. Ft. ]</b>	<b>9 Cubic Feet</b>
<b>Capacity [ Liters ]</b>	<b>255 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-0°C to -40°C</b>
<b>Temperature Range [ °F ]</b>	<b>+32°F to -40°F</b>

Dimensions (W x D x H)	
Exterior [ Inches ]	39" x 23.5" x 33.5"
Exterior [ Centimeters ]	99 x 59 x 85
Interior [ Inches ]	33.5" x 17.5" x 27.5"
Interior [ Centimeters ]	85 x 44 x 69

#### Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage	15 Amp Dedicated Line 2 Running Amps
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

#### Construction:

Temperature Control	Digital Display High / Low Alarm Min / Max Temperature
Alarm Relay	Dry Contacts
Defrost	Manual
Opening	Top Opening Lid
Lockable Lid/Door	No
Interior / Exterior	Steel
Insulation	Polyurethane Foam
Compressors	Hermetic Compressor
Primary Refrigerant(s)	R290
Air Circulation	Gravity Flow
Casters	4, Lockable
Shelving	None
Approximate BTU/Hour	600
Approximate Decibels	55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 1 Year Parts & 90 Days Labor