

MODEL DH-10RFC

LAB AND PHARMACY REFRIGERATOR FREEZER COMBO

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Total Capacity [Cu. Ft.] 10.2 Cu. Ft. (Total)
Capacity [Liters] 289 Liters (Total)

	REFRIGERATOR	FREEZER
Temperature Range [°C]	2°C to 8°C	-10°C to -26°C
Temperature Range [°F]	35°F to 46°F	14°F to -14.8°F

Dimensions	Width x Depth x Height	
Interior [Inches]	23" x 20" x 26"	19" x 18" x 19"
Interior [Centimeters]	60 x 51 x 67	50 x 46 x 51
Exterior [Inches]	25.5" x 25.75" x 72.75"	
Exterior [Centimeters]	75 x 65 x 185	


Compressor	Hermetic	Hermetic
Air Circulation	Direct	Forced Air
Defrost Method	Automatic	Manual
Porthole	One 1" Diameter	One 1" Diameter
Refrigerator Storage	3 Adjustable Shelves, Sliding Drawer	
Freezer Storage	3 Fixed Shelves, ABS Drawers	
Refrigerant	R600a	
Exterior	Painted Steel	
Interior	Stainless Steel	
Insulation	US EPA and SNAP approved	
Lockable Door	Yes, Lockable with Padlock	
Casters	4, with Leveling Feet	
Shipping Weight	300 lbs.	

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 1.85A
Power Consumption	3.11 kW/d
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Type	Microprocessor
Alarms	Two Digital Displays
Data Logging	High/Low Temp, Door Ajar
Min / Max Temperature	Power Failure, Low Battery
Alarm Relay	Sensor Failure, USB Failure
Backup Battery	Yes, USB Downloadable, PDF
	Yes, Display and Reset
	Yes, Dry Contacts
	Yes, Only Powers Alarm



Warranty 2 Year Parts & Labor