

MODEL DH-20RF

LAB AND PHARMACY REFRIGERATOR FREEZER COMBO

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

Total Capacity [Cu. Ft.] 18 Cu. Ft. (Total)
Capacity [Liters] 514 Liters (Total)

	REFRIGERATOR	FREEZER
Temperature Range [°C]	2°C to 8°C	0°C to -23°C
Temperature Range [°F]	35°F to 46°F	32°F to -10°F

Dimensions	Width x Depth x Height	
Interior [Inches]	24" x 28" x 26.5"	24" x 28" x 17.5"
Interior [Centimeters]	60 x 71 x 67	60 x 71 x 44
Exterior [Inches]	27.5" x 34" x 84"	
Exterior [Centimeters]	70 x 86 x 213	

Compressor	Hermetic	Hermetic
Air Circulation	Forced Air	Forced Air
Defrost Method	Automatic	Automatic
Porthole	1/2" Diameter	1/2" Diameter
Refrigerator Storage	2 Adjustable Shelves	
Freezer Storage	2 Adjustable Shelves	
Refrigerant	R290	
Exterior	Stainless Steel	
Interior	Stainless Steel	
Insulation	Polyurethane	
Lockable Door	Lockable with Padlock	
Casters	4, Locking	
Shipping Weight	300 lbs.	

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line	15 Amp Dedicated
Running Amperage	4.5A Refrigerator 7.5A Freezer
NEMA Configuration	NEMA 5-15
Plugs	1 Freezer, 1 Refrigerator
Certification(s)	UL Listed

Controller

Type	Microprocessor
Alarms	Two Digital Displays
Min / Max Temperature	High/Low Temp
Alarm Relay	Yes, Display and Reset
	Yes, Dry Contacts



Warranty 3 Year Parts & Labor.
 Additional 4 Years on compressor parts.