

## MODEL DH-20RFC

### LAB AND PHARMACY REFRIGERATOR FREEZER COMBO

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Total Capacity [ Cu. Ft. ] 18.33 Cu. Ft. (Total)**  
**Capacity [ Liters ] 519 Liters (Total)**

	REFRIGERATOR	FREEZER
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>	<b>-10°C to -26°C</b>
<b>Temperature Range [ °F ]</b>	<b>35°F to 46°F</b>	<b>14°F to -14.8°F</b>

Dimensions [ W x D x H ]		
Interior [ Inches ]	29.5" x 24"x 29"	25.5"x 22"x22.25"
Interior [ Centimeters ]	75 x 61 x 74	65 x 56 x 57
Exterior [ Inches ]	36" x 29.25" x 77.75"	
Exterior [ Centimeters ]	91 x 74 x 197	


Compressor	Hermetic	Hermetic
Air Circulation	Direct	Forced Air
Defrost Method	Automatic	Manual
Porthole	One 1" Diameter	One 1" Diameter
Refrigerator Storage	3 Adjustable Shelves, Sliding Drawer	
Freezer Storage	3 Fixed Shelves, ABS Drawers	
Refrigerant	R600a	
Exterior	Painted Steel	
Interior	Stainless Steel	
Insulation	US EPA and SNAP approved	
Lockable Door	Yes, with separate lock	
Casters	4, with Leveling Feet	
Shipping Weight	500 lbs.	

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   4.55A
Power Consumption	3.76 kW/d
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

#### Controller

Type	Microprocessor
Alarms	Two Digital Displays High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure
Data Logging	Yes, USB Downloadable, PDF
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts
Backup Battery	Yes, Only Powers Alarm



**Warranty** 2 Year Parts & Labor