

SO-LOW LAB & PHARMACY REFRIGERATOR/FREEZER MODEL DH-9RFDA

Refrigerator 4.6 cu. ft. (130 liters)-Freezer 4.4 cu. ft. (125 liters)

Upright – Temperature Range: 2°C/8°C (R) / -10°C/-30°C (F) Two Doors

Specifications:

* Capacity	9 cubic ft. / 255 liters
* Temperature Range	2°C to 8°C / -10°C to -30°C
* Dimensions	
* Exterior w x d x h	23.75" x 21.5" x 64" 60 x 54 x 162 (cm)
* Fridge Interior w x d x h	20" x 17.75" x 24.5"
* Freezer Interior w x d x h	20" x 17.75" x 24.5"

Features:

Temperature Control	Two - Digital display
Temperature Alarm	Hi/Low temperature alarm With audible, and visual alarm Alarm relay dry contacts
Doors	Two solid, locking
Legs	Four leveling feet
Shelves	Two adjustable / Refrigerator Two adjustable/ Freezer
Probe Port	3/8" access port hole, EA chamber
Stacking	Stacking Kit
Insulation	Urethane foam
Refrigerant:	R600A CFC free
Compressor	Two
Thermostat	Adjustable
Temperature	Temperature accuracy of $\pm 2-3^{\circ}\text{C}$
Refrigerator Defrost	Cycle, automatic
Freezer Defrost	Manual
Voltage: Domestic	115/60/1 – 1.2 amps (ref) 2.1 (frz)
Supply Plug	5-15P NEMA – 2 plugs
Shipping Weight	395 lbs.
Warranty-(Domestic)	Two year parts and labor

