

SO-LOW LAB & PHARMACY GLASS DOOR FREEZER MODEL DH20-23GD

23 Cubic Ft. Upright – Temperature Range: -15°C to -25°C

Glass Door

Specifications:

* Capacity	23 cubic ft.
* Temperature Range	-15°C to -25°C
* Dimensions	
* Exterior w x d x h	27.25" x 34" x 81" 69 x 86 x 206 (cm)
* Interior w x d x h	24" x 28.5" x 55" 61 x 72 x 140 (cm)

Features:

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Doors	Single glass, locking
Interior finish	White powder coat
Casters	Four, Locking
Shelves	Five, Adjustable
Interior Lighting	One fluorescent
Insulation	Polyurethane foam
Refrigerant:	R404A CFC free
Compressor	1/2 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	10.9/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Two year parts and labor Compressor additional 3 years, no labor.

