

SO-LOW LAB & PHARMACY GLASS DOOR FREEZER MODEL DH20-27GD

27 Cubic Ft. Upright – Temperature Range: -15°C to -25°C

Glass Door

Specifications:

* Capacity	27 cubic ft.
* Temperature Range	-15°C to -25°C
* Dimensions	
* Exterior w x d x h	30" x 34" x 81" 76 x 86 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

Features:

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Doors	Single glass, locking
Interior finish	White powder coat, s/s floor
Casters	Four
Shelves	Five
Interior Lighting	One 20-watt fluorescent
Insulation	Urethane foam
Refrigerant:	R404A CFC free
Compressor	3/4 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	13.8/20
Supply Plug	5-20P NEMA UL listed
Shipping Weight	450 lbs.
Warranty-(Domestic)	Two year parts and labor Compressor additional 3 years, no labor.

