

SO-LOW LAB & PHARMACY GLASS DOOR FREEZER MODEL DH20-72GD

72 Cubic Ft. Upright – Temperature Range: -20°C to -26°C

Glass Doors

Specifications:

* Capacity	72 cubic ft.
* Temperature Range	-20°C to -26°C
* Dimensions	
* Exterior w x d x h	75" x 34" x 81" 132 x 86 x 206 (cm)
* Interior w x d x h	72" x 28.5" x 55" 183 x 72 x 140 (cm)

Features:

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Doors	Three glass, locking
Interior finish	White powder coat, s/s floor
Casters	Four
Shelves	Fifteen
Interior Lighting	One 20-watt fluorescent
Insulation	Urethane foam
Refrigerant:	R404A CFC free
Compressor	3/4 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, automatic
Voltage: Domestic	230/115/60/1
Running Amps/Breaker	15/20
Supply Plug	L14-20 NEMA UL listed
Shipping Weight	800 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.

