

MODEL DH25-23SD-SS

LAB AND PHARMACY FREEZER -18°C TO -25°C

Capacity [Cu. Ft.]	23.32 Cu. Ft.
Capacity [Liters]	660 Liters
Temperature Range [°C]	-18°C to -25°C
Temperature Range [°F]	0°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	27.25" x 33.75" x 81.25"
Interior [Inches]	24.25" x 28.5" x 61.75"
Exterior [Centimeters]	69 x 85 x 206
Interior [Centimeters]	62 x 72 x 157

Compressor	3/4 HP
Air Circulation	Forced Air
Defrost Method	Adaptive defrost
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Stainless Steel
Interior Light	LED
Porthole	1/2" Diameter (LH Side)
Casters	4, Front Locking
Storage	3 Adjustable Shelves
Shipping Weight	450 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6.34A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

Warranty 3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

