

MODEL DH25-23SD

LAB AND PHARMACY FREEZER -18°C TO -25°C

Capacity [Cu. Ft.]	23.32 Cu. Ft.
Capacity [Liters]	660 Liters
Temperature Range [°C]	-18°C to -25°C
Temperature Range [°F]	0°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	27.25" x 33.5" x 81"
Interior [Inches]	24.25" x 28.5" x 61.75"
Exterior [Centimeters]	69 x 85 x 206
Interior [Centimeters]	61 x 72 x 157

Compressor	3/4 HP
Air Circulation	Forced Air
Defrost Method	Auto Defrost
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Interior Light	LED
Porthole	1/2" Diameter (LH Side)
Casters	4, Front Locking
Storage	5 Adjustable Shelves
Shipping Weight	450 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6.34A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.