

## MODEL DH25-27SD-SS

### LAB AND PHARMACY FREEZER -18°C TO -25°C

|                                 |                       |
|---------------------------------|-----------------------|
| <b>Capacity [ Cu. Ft. ]</b>     | <b>26 Cu. Ft.</b>     |
| <b>Capacity [ Liters ]</b>      | <b>736 Liters</b>     |
| <b>Temperature Range [ °C ]</b> | <b>-18°C to -25°C</b> |
| <b>Temperature Range [ °F ]</b> | <b>0°F to -13°F</b>   |

|                          |                        |
|--------------------------|------------------------|
| Dimensions               | Width x Depth x Height |
| Exterior [ Inches ]      | 30" x 33.75" x 81.25"  |
| Interior [ Inches ]      | 27" x 28.5" x 61.75"   |
| Exterior [ Centimeters ] | 76 x 85 x 206          |
| Interior [ Centimeters ] | 69 x 72 x 157          |

|                 |                         |
|-----------------|-------------------------|
| Compressor      | 3/4 HP                  |
| Air Circulation | Forced Air              |
| Defrost Method  | Adaptive defrost        |
| Refrigerant     | R290                    |
| Exterior        | Stainless Steel         |
| Interior        | Stainless Steel         |
| Insulation      | Polyurethane Foam       |
| Lockable Door   | 1 Door, Keyed Lock      |
| Door Material   | Stainless Steel         |
| Interior Light  | LED                     |
| Porthole        | 1/2" Diameter (LH Side) |
| Casters         | 4, Front Locking        |
| Storage         | 3 Adjustable Shelves    |
| Shipping Weight | 550 lbs.                |

#### Electrical

|                         |                            |
|-------------------------|----------------------------|
| Voltage:                | 115 Volts 60 Hertz 1 Phase |
| Amperage Line   Running | 15 Amp Dedicated   10A     |
| NEMA Configuration      | NEMA 5-15                  |
| Certification(s)        | UL Listed                  |

#### Controller

|                       |                         |
|-----------------------|-------------------------|
| Display Type          | Microprocessor, Digital |
| Alarms                | High/Low Temp           |
| Min / Max Temperature | Display and Reset       |
| Alarm Relay           | Dry Contacts            |

#### Warranty

3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

