

## MODEL DH25-49SD-SS

### LAB AND PHARMACY FREEZER -15°C TO -25°C

<b>Capacity [ Cu. Ft. ]</b>	<b>46 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>1387 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-15°C to -25°C</b>
<b>Temperature Range [ °F ]</b>	<b>5°F to -13°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	52" x 33 3/4" x 81"
Interior [ Inches ]	49" x 26" x 61 5/8"
Exterior [ Centimeters ]	132 x 81 x 206
Interior [ Centimeters ]	120 x 72 x 156

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Adaptive defrost
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Doors	Two Doors, Keyed Lock
Door Material	Insulated Steel
Interior Light	LED
Porthole	1/2" Diameter Left Side
Casters	4, Locking
Storage	6 Adjustable Shelves
Shipping Weight	670 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	20 Amp Dedicated   11A
NEMA Configuration	NEMA 5-20
Certification(s)	UL Listed

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Alarm Relay	Dry Contacts

#### Warranty

3 Years Parts & Labor.  
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

