

SO-LOW LAB & PHARMACY FREEZER *MODEL DH29-23SD*

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

24 Cubic Ft. (679 liters) Solid Door Upright

Temperature Range: -18°C to -25°C / 0°F to -13°F

Specifications:

* Capacity	24 cubic ft. / 679 liters
* Temperature Range	-18°C to -25°C / 0°F to -13° F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	26" x 34" x 81.5" 66 x 86 x 207 (cm)
* Interior w x d x h	22" x 28" x 60" 56 x 71 x 152 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Door	Single solid – lockable
Casters	heavy duty, swivel
Shelves	3 epoxy coated adjustable
Interior Lights	One
Insulation	Polyurethane CFC Free
Refrigerant:	R-290
Air Circulation	Forced air
Refrigerant Control	Expansion valve
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Timed automatic
Voltage:	115/60/1
Running Amps/Breaker	7.1/15
Supply Plug	5-15 NEMA
Shipping Weight	425 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

