

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## SO-LOW LAB & PHARMACY FREEZER

### MODEL DH29-34SD

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

**34 Cubic Ft. (963 liters) Upright**

**Temperature Range: -18°C to -25°C / 0°F to -25°F**  
**Solid Door**

#### Specifications:

\* Capacity 34 cubic ft. / 963 liters  
\* Temperature Range -18°C to -25°C / 0°F to -13°F  
**\*When ordering please specify Celsius or Fahrenheit Setting\***

#### \* Dimensions

\* Exterior w x d x h 35" x 34" x 81.5"  
89 x 86 x 207 (cm)

\* Interior w x d x h 31" x 28" x 60"  
79 x 71 x 152 (cm)

#### Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Door	Single solid – lockable
Casters	heavy duty, swivel
Shelves	3 epoxy coated adjustable
Interior Lights	One 40W
Insulation	Polyurethane CFC Free
Refrigerant:	R-404A
Air Circulation	Forced air
Refrigerant Control	Expansion valve
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Timed automatic
Voltage:	115/60/1
Running Amps/Breaker	7.8/15
Supply Plug	5-15 NEMA
Shipping Weight	450 lbs.
Warranty-(Domestic)	Three year parts and labor Compressor additional 2 years, no labor.



Specifications subject to change without notice.