

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB & PHARMACY FREEZER

MODEL DH29-46SD

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

49 Cubic Ft. (1387 liters) Upright

Temperature Range: -18°C to -25°C / 0°F to -13°F

2 Solid Doors

Specifications:

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	-18°C to -25°C / 0°F to -13° F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	52" x 34" x 81.5" 132 x 86 x 207 (cm)
* Interior w x d x h	48" x 28" x 60" 122 x 71 x 152 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Doors	Two doors solid - lockable
Casters	Heavy duty, swivel
Shelves	6 epoxy coated adjustable
Interior Lights	Two 40W
Insulation	Polyurethane CFC Free
Refrigerant:	R-404A
Compressor	3/4 hp
Air Circulation	Forced air
Refrigerant Control	Expansion valve
Evaporator	Fin and tube, forced air
Defrost	Timed automatic
Voltage:	115/60/1
Running Amps/Breaker	12/20
Supply Plug	5-20 NEMA
Shipping Weight	650 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice.