

## *SO-LOW LABORATORY FREEZER MODEL DH30-74SDSS*

**72 Cubic Ft. (2095 liters) Upright –Temperature Range: 0°C to -20°C / 32°F to -4°F  
Three Hinged Solid Doors**

### *Specifications:*

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	0°C to -20°C / 32°F to -4° F
<b>*When ordering please specify Celsius or Fahrenheit Setting*</b>	
* Dimensions	
* Exterior w x d x h	75" x 33.75" x 81" 190 x 85 x 206 (cm)
* Interior w x d x h	72" x 27" x 55" 183 x 69 x 140 (cm)

### *Features:*

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Exterior	Stainless Steel
Interior	Stainless Steel
Doors	Three solid – hinged – key lock Positive seal, self closing
Shelves	Nine adjustable
Casters	Four, Locking
Interior Lighting	Two
Insulation	Polyurethane CFC Free
Refrigerant	R290 CFC free
Compressor	3/4 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, automatic
Voltage: Domestic	208-230
Running Amps/Breaker	13.8/20
Supply Plug	NEMA 14-20P UL listed
Shipping Weight	900 lbs.
Warranty-(Domestic)	Three year parts and labor Compressor additional 2 years, no labor.

