

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-10GD

**10 Cubic Ft. Upright – Temperature Range: 2°C to 8°C
Glass Door**

Specifications:

* Capacity	10 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	24" x 26" x 65" 60 x 66 x 165 (cm)
* Interior w x d x h	20.75" x 21" x 37.75" 52 x 53 x 96 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Door	Single glass, locking
Access Port	1/2" Diameter Top Back Wall
Interior finish	White coated steel
Casters	Four locking
Shelves	Three Adjustable
Interior Lighting	One fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/5 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	300 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

