

## *SO-LOW LABORATORY REFRIGERATOR MODEL DH4-17SGD*

*"SETTING TEMPERATURE STANDARDS IN EVERY LAB"*

**17 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**Two Glass Sliding Doors**

### *Specifications:*

* Capacity	17 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	36" x 24" x 71" 91 x 60 x 180 (cm)
* Interior w x d x h	33 x 21" x 46" 83 x 53 x 116 (cm)

### *Features:*

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Exterior	White epoxy powder-coated
Interior	White epoxy powder-coated
Doors	Double glass – sliding Positive seal, self closing
Casters	Four, locking
Shelves	Four
Interior Lighting	One fluorescent light
Insulation	Foamed in place polyurethane
Refrigerant:	R134A
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.8/15
Supply Plug	5-15P NEMA
Shipping Weight	450 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

