

## MODEL DH4-23GD-SS

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>23 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>651 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	27.25" x 33.75" x 81"
Interior [ Inches ]	24" x 28.5" x 55"
Exterior [ Centimeters ]	69 x 85 x 206
Interior [ Centimeters ]	61 x 72 x 140

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	3 Adjustable Shelves
Shipping Weight	440 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   3.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 7 Years Parts, Labor & Compressor Part Warranty.