

## *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-23SD*

**23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Solid Door**

### **Specifications:**

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	27" x 34" x 81" 68 x 86 x 206 (cm)
* Interior w x d x h	24" x 28.5" x 55" 61 x 72 x 140 (cm)

### **Features:**

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter Top Back Wall
Doors	Single Solid, locking
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four, lockable
Shelves	Five, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice