

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-27GD

**27 Cubic Ft. (764 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Glass Door**

Specifications:

* Capacity	27 cubic ft. / 764 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	30" x 32" x 81" 76 x 81 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Door	Single glass, locking
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four, locking
Shelves	Five, Adjustable
Interior Lighting	One
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.0 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

