

MODEL DH4-27GD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 27 Cu. Ft.
Capacity [Liters] 764 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions Width x Depth x Height
Exterior [Inches] 30" x 32" x 81"
Interior [Inches] 27" x 28.5" x 55"
Exterior [Centimeters] 76 x 81 x 206
Interior [Centimeters] 69 x 72 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Door 1 Door, Keyed Lock
Door Material Glass
Interior Light Adjustable LED
Porthole 1/2" Diameter (Top Back Wall)
Casters 4, Locking
Storage 5 Adjustable Shelves
Shipping Weight 440 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 3.0A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed
Energy Star Rated

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.