

## *SO-LOW CHROMATOGRAPHY REFRIGERATOR MODEL DH4-27GDCH*

**27 Cubic Ft. (764 liters) Upright – Temperature Range: 2°C to 8°C**

### **Glass Door**

#### *Specifications:*

* Capacity	27 cubic ft. / 764 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	30" x 32" x 81" 76 x 81 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

#### *Features:*

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Interior	White powder-coated steel
Interior Outlet	115V Duplex 5A
Access Ports	LH & RH Ends of cabinet 2" diameter
Door	Single glass – hinged Positive seal, self closing
Shelves	Five
Casters	Four, Locking
Interior Lighting	One fluorescent light
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps	5.8/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

