

MODEL DH4-27SD-SS

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	27 Cu. Ft.
Capacity [Liters]	764 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F
Dimensions	Width x Depth x Height
Exterior [Inches]	30" x 33.75" x 81"
Interior [Inches]	27" x 28.5" x 55"
Exterior [Centimeters]	76 x 85 x 206
Interior [Centimeters]	69 x 72 x 140
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Stainless Steel
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	3 Adjustable Shelves
Shipping Weight	440 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 3.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.