

SO-LOW CHROMATOGRAPHY REFRIGERATOR MODEL DH4-38SGDCH

38 Cubic Ft. (1076 liters) Upright – Temperature Range: 2°C to 8°C

Two Glass Sliding Doors

Specifications:

* Capacity	38 cubic ft. / 1076 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	43.5" x 32" x 81" 110 x 81 x 198 (cm)
* Interior w x d x h	40.5 x 27" x 55" 102 x 68 x 140 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Interior	White powder-coated steel
Interior Outlet	115V Duplex 5A
Ports	Access ports on LH & RH Ends of cabinet – 2" diameter
Doors	Double glass – sliding Positive seal, self closing Triple pane thermal glass
Shelves	Ten , epoxy coated wire
Interior Lighting	One fluorescent light
Insulation	Polyurethane foam
Refrigerant:	R134A
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps	8.8/15 Amp circuit required
Supply Plug	5-15P NEMA
Shipping Weight	525 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

