

MODEL DH4-49GD-SS

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	49 Cu. Ft.
Capacity [Liters]	1387 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F

Dimensions	Width x Depth x Height
Exterior [Inches]	52" x 33.75" x 81"
Interior [Inches]	47.5" x 28.5" x 55"
Exterior [Centimeters]	132 x 85 x 206
Interior [Centimeters]	120 x 72 x 140

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Door	2 Doors, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	6 Adjustable Shelves
Shipping Weight	610 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 5.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor.
 Additional 4 Years on Compressor Parts, No Labor.