

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-49GD

49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Two Swinging Glass Doors

Specifications:

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	52" x 32" x 81" 132 x 81 x 206 (cm)
* Interior w x d x h	49" x 28.5" x 55" 125 x 72 x 140 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter Top Back Wall
Doors	Two hinged glass, locking
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four, Locking
Shelves	Ten, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	610 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.

