

MODEL DH4-49SD-SS

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 49 Cu. Ft.
Capacity [Liters] 1387 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions Width x Depth x Height
Exterior [Inches] 52" x 33.75" x 81"
Interior [Inches] 47.5" x 28.5" x 55"
Exterior [Centimeters] 132 x 85 x 206
Interior [Centimeters] 120 x 72 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Stainless Steel
Interior Stainless Steel
Insulation Polyurethane Foam
Lockable Door 2 Doors, Keyed Lock
Door Material Insulated Stainless Steel
Interior Light Adjustable LED
Porthole 1/2" Diameter
Casters 4, Locking
Storage 6 Adjustable Shelves
Shipping Weight 610 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 5.0A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed
Energy Star Rated

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.