

MODEL DH4-49SD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 49 Cu. Ft.
Capacity [Liters] 1387 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions Width x Depth x Height
Exterior [Inches] 52" x 32" x 81"
Interior [Inches] 47.5" x 28.5" x 55"
Exterior [Centimeters] 132 x 81 x 206
Interior [Centimeters] 120 x 72 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Doors Two Doors, Keyed Lock
Door Material Insulated Steel
Interior Light Adjustable LED
Porthole 1/2" Diameter (Top Back Wall)
Casters 4, Locking
Storage 10 Adjustable Shelves
Shipping Weight 610 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 5.0A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed
Energy Star Rated

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.