

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-49SDSS

**49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Two Solid Doors**

Specifications:

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	52" x 33.75" x 81"
	132 x 85 x 206 (cm)
* Interior w x d x h	49" x 28.5" x 55"
	125 x 72 x 140 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter Top Back Wall
Doors	Two solid locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four, Locking
Shelves	Six, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	610 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.

