

## MODEL DH4-72GD-SS

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

**Capacity [ Cu. Ft. ]** 72 Cu. Ft.  
**Capacity [ Liters ]** 2095 Liters  
**Temperature Range [ °C ]** 2°C to 8°C  
**Temperature Range [ °F ]** 36°F to 46°F

**Dimensions**  
Exterior [ Inches ] 75" x 33.75" x 81"  
Interior [ Inches ] 72" x 27" x 55"  
Exterior [ Centimeters ] 190 x 85 x 206  
Interior [ Centimeters ] 183 x 69 x 140

**Compressor** Hermetic Compressor  
**Air Circulation** Forced Air  
**Defrost Method** Cycle  
**Refrigerant** R290  
**Exterior** Stainless Steel  
**Interior** Stainless Steel  
**Insulation** Polyurethane Foam  
**Lockable Door** 3 Doors, Keyed Lock  
**Door Material** Glass  
**Interior Light** Adjustable LED  
**Porthole** 1/2" Diameter  
**Casters** 4, Locking  
**Storage** 9 Adjustable Shelves  
**Shipping Weight** 850 lbs.

#### Electrical

**Voltage:** 115 Volts 60 Hertz 1 Phase  
**Amperage Line | Running** 15 Amp Dedicated | 7.2A  
**NEMA Configuration** NEMA 5-15  
**Certification(s)** UL Listed  
Energy Star Rated

#### Controller

**Display Type** Microprocessor, Digital  
**Alarms** High/Low Temp  
**Min / Max Temperature** Display and Reset  
**Alarm Relay** Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 7 Years Parts, Labor & Compressor Part Warranty.