

SO-LOW CHROMATOGRAPHY REFRIGERATOR MODEL DH4-74GDCH-SS

**72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C
Three Glass Hinged Doors**

Specifications:

| | |
|----------------------|---|
| * Capacity | 72 cubic ft. / 2095 liters |
| * Temperature Range | 2°C to 8°C |
| * Dimensions | |
| * Exterior w x d x h | 75" x 33.75" x 81" 190 x 85 x 206 (cm) |
| * Interior w x d x h | 72" x 27" x 55" 183 x 69 x 140 (cm) |

Features:

| | |
|-------------------------------|--|
| Temperature Control | Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts |
| Interior Outlet Ports | 115V Duplex 5A Access ports on LH & RH Ends of cabinet – 2" diameter |
| Exterior Interior Doors | Stainless Steel Stainless Steel Three glass – hinged – key lock Positive seal, self closing |
| Shelves | Nine, Adjustable |
| Casters | Four, Locking |
| Interior Lighting | Two fluorescent lights |
| Insulation | Polyurethane CFC Free |
| Refrigerant | R134A CFC free |
| Compressor | 1/2 hp |
| Evaporator | Fin and tube, forced air |
| Defrost | Cycle, Auto |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 9.9/15 |
| Supply Plug | NEMA 5-15P UL listed |
| Shipping Weight | 850 lbs. |
| Warranty-(Domestic) | Two year parts and labor Compressor additional 4 years, no labor. |

