

## SO-LOW CHROMATOGRAPHY REFRIGERATOR

### MODEL DH4-74GDCH

**72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C  
Three Glass Hinged Doors**

#### Specifications:

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	75" x 32" x 81" 190 x 81 x 206 (cm)
* Interior w x d x h	72" x 27" x 55" 183 x 69 x 140 (cm)

#### Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Interior	White epoxy powder-coated steel
Interior Outlet Ports	115V Duplex 5A Access ports on LH & RH Ends of cabinet – 2" diameter
Doors	Three glass – hinged w/key lock Positive seal, self closing Triple pane thermal glass
Shelves	Fifteen, Adjustable
Interior Lighting	Two fluorescent
Insulation	Polyurethane
Refrigerant:	R134A CFC free
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps	9.9/15 Amp circuit required
Supply Plug	5-15 NEMA UL listed
Shipping Weight	750 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.

