

## *SO-LOW LABORATORY REFRIGERATOR MODEL DH4-74GDSS*

**72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C  
Three Glass Hinged Doors**

### *Specifications:*

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	75" x 33.75" x 81" 190 x 85 x 206 (cm)
* Interior w x d x h	72" x 27" x 55" 183 x 69 x 140 (cm)

### *Features:*

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Exterior	Stainless Steel
Interior	Stainless Steel
Doors	Three glass – hinged – key lock Positive seal, self closing
Shelves	Nine, Adjustable
Casters	Four, Locking
Interior Lighting	Two fluorescent lights
Insulation	Polyurethane CFC Free
Refrigerant	R134A CFC free
Compressor	1/2 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	9.9/15
Supply Plug	NEMA 5-15P UL listed
Shipping Weight	850 lbs.
Warranty-(Domestic)	Two year parts and labor Compressor additional 3 years, no labor.

