

MODEL DH4-74SD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 72 Cu. Ft.
Capacity [Liters] 2095 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions Width x Depth x Height
Exterior [Inches] 75" x 32" x 81"
Interior [Inches] 72" x 27" x 55"
Exterior [Centimeters] 190 x 81 x 206
Interior [Centimeters] 183 x 69 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Doors Three Doors, Keyed Lock
Door Material Insulated Steel
Interior Light Adjustable LED
Porthole 1/2" Diameter (Top Back Wall)
Casters 4, Locking
Storage 15 Adjustable Shelves
Shipping Weight 720 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 7.2A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed
Energy Star Rated

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Yes, Display and Reset
Alarm Relay Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.