

MODEL DH4-80PT

PASS-THRU LAB REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	80 Cu. Ft.
Capacity [Liters]	2265 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F
Dimensions	Width x Depth x Height
Exterior [Inches]	78" x 37.5" x 81.5"
Interior [Inches]	75" x 30" x 59"
Exterior [Centimeters]	198 x 95 x 207
Interior [Centimeters]	190 x 76 x 149
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Door	6 Doors, Keyed Lock
Door Material	Glass (Front & Back Doors) Option for Solid Door(s)
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	9 Adjustable Shelves
Shipping Weight	850 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 14.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

Warranty 3 Years Parts & Labor.
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

