

MODEL DHH-30RFFMS

FLAMMABLE MATERIAL STORAGE REFRIGERATOR / FREEZER COMBINATION UNIT

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

Total Capacity [Cu. Ft.] 30 Cu. Ft. (Total)
Capacity [Liters] 849 Liters (Total)

	REFRIGERATOR	FREEZER
Temperature Range [°C]	1°C to 10°C	-15°C to -25°C
Temperature Range [°F]	33°F to 55°F	5°F to -13°F

Dimensions	Width x Depth x Height	
Interior [Inches]	29" x 22.25" x 12"	32" x 22.5" x 50.25"
Interior [Centimeters]	73 x 56 x 30	81 x 57 x 127
Exterior [Inches]	35.75" x 35.5" x 72.5"	
Exterior [Centimeters]	90 x 90 x 184	

Compressor	1 Each, Hermetic Compressor
Refrigerator Air Circulation	Gravity Flow
Refrigerator Defrost Method	Manual
Refrigerator Storage	3 Adjustable Shelves
Freezer Air Circulation	Gravity Flow
Freezer Defrost Method	Manual
Freezer Storage	N/A
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Leveling Legs	4
Shipping Weight	400 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Mechanical Thermostat
--------------	-----------------------



» Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.

Warranty 1 Year Parts & Labor