

SO-LOW EXPLOSION PROOF MANUAL DEFROST REFRIGERATOR/FREEZER MODEL DHH-30RFX

Meets NFPA and OSHA guidelines for Hazardous Locations, Class 1, Division II, Group C & D.

Refrigerator 25 cu. ft. (707 liters)-Freezer 5.0 cu. ft. (141 liters)

Upright – Temperature Range: 1°C to 10°C(R) / -18°C to -25°C (F)

Specifications:

* Capacity	30 cubic ft. / 849 liters
* Temperature	1°C to 10°C / -18°C to -25°C
* Dimensions	
* Exterior w x d x h	35 ¾" x 36 ¼" x 77 ½" 90.81 x 92.08 x 196.85 (cm)
* Interior w x d x h	29" x 22.25" x 12" - Freezer 73.66 x 56.52 x 30.48 (cm) 32" x 22.5" x 50.25 - Refrigerator 81.28 x 57.15 x 127.64 (cm)

Features:

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	3 adjustable – Fridge only
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-404A
*Compressor	1/3 hp
*Air Circulation	Gravity
*Evaporator	Hot wall
*Defrost	Manual
*Voltage	120 VAC
*Running Amps	6
*Supply Plug	None / Unit must be hardwired
*Shipping weight	408 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.

