

## MODEL DHH20-14MDP

### MANUAL DEFROST FREEZER -15°C TO -25°C

<b>Capacity [ Cu. Ft. ]</b>	<b>14 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>396 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-15°C to -25°C</b>
<b>Temperature Range [ °F ]</b>	<b>5°F to -13°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	23.63" x 24.38" x 73"
Interior [ Inches ]	19" x 17.5" x 63.75"
Exterior [ Centimeters ]	60 x 62 x 185
Interior [ Centimeters ]	48 x 44 x 161

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R290
Exterior	Powder Coating
Interior	Molded ABS
Insulation	Polyurethane Foam
Lockable Door	Keyed
Porthole	3/4" Diameter
Leveling Legs	4
Storage	6 Fixed Shelves
Shipping Weight	225 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   3.5A
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Alarm Relay	Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor