

MODEL DHH20-14SDFMS

FLAMMABLE MATERIAL STORAGE FREEZER -15°C TO -25°C

Capacity [Cu. Ft.] 14 Cu. Ft.
Capacity [Liters] 396 Liters
Temperature Range [°C] -15°C to -25°C
Temperature Range [°F] 5°F to -13°F

Dimensions Width x Depth x Height
 Exterior [Inches] 23.75" x 24.75" x 73"
 Interior [Inches] 19" x 17.5" x 63.75"
 Exterior [Centimeters] 60 x 63 x 185
 Interior [Centimeters] 48 x 44 x 162

Compressor Hermetic Compressor
 Air Circulation Gravity Flow
 Defrost Method Manual
 Refrigerant R290
 Exterior Painted Steel
 Interior Painted Steel
 Insulation High Density Urethane
 Lockable Door 1 Door, Keyed Lock
 Door Material Insulated Steel
 Porthole 3/4" Diameter
 Leveling Legs 4
 Storage 6 Fixed Shelves
 Shipping Weight 195 lbs.

Electrical
 Voltage: 115 Volts 60 Hertz 1 Phase
 Amperage Line | Running 15 Amp Dedicated | 3.5A
 NEMA Configuration NEMA 5-15
 Certification(s) ETL Listed
 Energy Star Rated

Controller
 Display Type Microprocessor, Digital
 Alarms High/Low Temp
 Min / Max Temperature Display and Reset
 Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.**

Warranty 2 Years Parts & Labor.

Specifications subject to change without notice